



ZirChrom®

# A Simple and Sensitive HPLC Method for the Detection and Quantitation Of STI-571 And Its Main Metabolite N-Desmethyl-STI

Clayton McNeff, Ph.D. and Steven Rupp  
ZirChrom Separations, Inc.

## Technical Bulletin # 302

An isocratic online enrichment HPLC-assay was developed allowing for the simple and fast separation and quantitation of STI-571 and its main metabolite N-Desmethyl-STI in plasma, urine, cerebrospinal fluid, culture media and cell preparations in various concentrations using UV-detection at 260 nm. The analytical procedure consists of an online concentration of STI-571 and N-Desmethyl-STI in the HPLC system followed by the elution on a ZirChrom®-PBD analytical column.




Figure 1. Structures of STI-571 and N-Desmethyl-STI

### Introduction

STI-571 (Imatinib mesylate, Glivec™), a 4-[(4-methyl-1-piperazinyl)methyl]-N-[4-[1-3]methyl-3-[[4-(3-pyridinyl)-2-pyrimidinyl]amino]-phenyl]benzamide methane-sulfonate derivative, acts as an inhibitor of the abl tyrosine kinase, platelet derived growth factor receptor (PDGFR), stem cell factor receptor (c-kit, steel factor receptor, CD117) and ARG tyrosine kinases. The specific blockade of the bcr-abl oncoprotein has been associated with significant antileukemic activity in patients with chronic-myeloid-leukemia (CML) and Philadelphia-positive-acute-lymphatic-leukemia (Ph+ALL).

This analytical method was developed by Universitätsklinikum Carl Gustav Carus an der Technischen Universität (Dresden, Germany) to perform pharmacokinetic measurements of STI-571 and N-Desmethyl-STI in patient samples (plasma, urine, cerebrospinal fluid) and for kinetic measurements of intracellular STI-571 and N-Desmethyl-STI following in-vitro incubation [1]. This method utilizes UV detection but may also be adapted to electrochemical detection to enable lower detection limits.

### Experimental

A mixture of STI-571 and N-Desmethyl-STI was separated at room temperature using a ZirChrom®-PBD guard column, a ZirChrom®-PBD analytical column and UV detection. This analytical method includes an online-enrichment system incorporating another ZirChrom®-PBD guard column and an electric motor driven switching valve (refer to [1] for schematic valve switching graph). The separation conditions were as follows:

#### Enrichment

Guard Column: ZirChrom®-PBD, 10 mm x 4.6 mm i.d. Guard Insert (Part Number: ZR03-G40; set of 3); Guard Insert Holder (Part Number 850-00)  
Mobile Phase: 45/35/20 (v/v) 0.1 M dibasic potassium phosphate/water/methanol  
Temperature: Uncontrolled

#### Analytical

Guard Column: ZirChrom®-PBD, 10 mm x 4.6 mm i.d. Guard Insert (Part Number: ZR03-G40; set of 3); Guard Insert Holder (Part Number 850-00)  
Column: ZirChrom®-PBD, 50 mm x 4.6 mm i.d. (Part Number ZR03-0546)  
Mobile Phase: 60/40 (v/v) 10 mM dibasic potassium phosphate, 90 mM monobasic potassium phosphate/methanol  
Temperature: Uncontrolled  
Flow Rate: 0.4 ml/min.  
Injection: 50 µl  
Detection: UV at 260 nm

Time of analysis is 40 minutes including the enrichment time of 5 minutes. The UV detection limit is 10 ng/ml in plasma, CSF (cerebrospinal fluid), culture medium (RPMI) and 25 ng/ml in urine for both STI-571 and N-Desmethyl-STI.




Figure 2: Chromatogram of patient plasma under STI-571 treatment with 5308 ng/ml STI-571 (RT ≈ 10.0 min.) and 988 ng/ml N-Desmethyl-STI (RT ≈ 30.0 min.).

This method can be tailored to your specific application needs. ZirChrom method developers can help to optimize and transfer this method to your site. Please contact ZirChrom technical support at 1-866-STABLE-1 or [support@zirchrom.com](mailto:support@zirchrom.com) for details.

#### Acknowledgements

Universitätsklinikum Carl Gustav Carus an der Technischen Universität (Dresden, Germany)

#### References

- [1] E. Schleyer, S. Pursche, C.H. Köhne, U. Schuler, U. Renner, H. Gschaidmeier, J. Freiberg-Richter, T. Leopold, A. Jenke, M. Bonin, T. Bergemann, P. le Coutre, M. Gruner, M. Bornhäuser, O.G. Ottmann, G. Ehninger, J. Chromatogr. B, 799, (2004) 23-36.

**ZirChrom Separations, Inc.**  
617 Pierce Street, Anoka, MN 55303  
1-866-STABLE-1  
[support@zirchrom.com](mailto:support@zirchrom.com)

Visit [www.zirchrom.com](http://www.zirchrom.com) for more application notes using ultra-stable, high efficiency ZirChrom columns.