

The Preparation of Buffers with **ZirChrom**'s Buffer Wizard

<http://www.zirchrom.com>

What is pH ?

- A scale for measuring solution acidity/basicity

$$\text{pH} = -\log a_{\text{H}^+} \cong -\log [\text{H}^+]$$

- In dilute solution, H^+ activity is essentially equal to the H^+ ion concentration.
- pH in water is commonly between 0 and 14

What is a pH buffer?

- A solution that **resists** changes in pH when **small** amounts of acids or bases are added to the solution or when the solution is diluted.
- The buffer solution is commonly prepared from a weak acid and its conjugate base or from a weak base and its conjugate acid.

How does a pH buffer resist changes in pH?

- A pH buffer contains species that can neutralize small amounts of acids & bases.

pH and Liquid Chromatography

- pH of the eluent may be changed by the sample being analyzed.
- The pH changes affect acid-base equilibria in the system.
- Shifts in the acid-base equilibria may affect reproducibility, selectivity, and peak shape for ionizable compound.
- High & low pH can also destabilize silica-based columns.

Why Control/Adjust pH?

- To improve (increase/decrease) retention.
- To adjust band spacing.
- To improve peak shape (tailing).
- To make results reproducible.
- To improve detection (UV, MS).

Buffer Capacity: A Key Concept

- Ability of a buffer to resist changes in pH upon addition of an acid or a base.
- Defined as the number of moles of strong acid or base per liter required to add to the solution to produce a unit change in pH.

$$\beta = dC_b / dpH = - dC_a / dpH$$

Acetic Acid

A Simple Monoprotic Acid-Base Equilibrium




$$[H^+] = K_a \cdot \frac{[HAc]}{[Ac^-]} \quad pH = pK_a + \log \frac{[Ac^-]}{[HAc]}$$


Henderson - Hasselbalch Eq.

$$pK_a = 4.76$$

Titration of 20 mL 0.01 M HAc with 0.01 M NaOH



Buffer Capacity of 0.01 M HAc



Phosphoric Acid

A Triprotic Acid-Base System




$$pK_{a,1} = 2.15$$


$$pK_{a,2} = 7.20$$

$$pK_{a,3} = 12.15$$

Titration of 20 mL 0.01 M H_3PO_4 with 0.01 M NaOH




Common pH Buffers



Acetic Acid + Ammonium Hydroxide


A Mixture of Two Monoprotic Acid-Base Systems




$$pK_{a,HAc} = 4.76$$

$$pK_{a,NH_4OH} = 9.24$$

Titrating 20 mL 0.01 M HAc + 0.01 M NH₄OH
with 0.01 M NaOH



Buffer Capacity of 0.01 M HAc + 0.01 M NH₄OH




Buffer Action upon Dilution

- Weak acid + Weak base
0.085 M HAc + 0.015 NaAc
Initial pH: 4.0
pH after 10 times dilution: 4.0
- Strong acid
0.0001 M HCl
Initial pH: 4.0
pH after 10 times dilution: 5.0


Key Issues in Buffer Preparation for HPLC

- Adequate buffer capacity—pick buffer with pK_a closest to the desired pH.
- Buffer solubility & compatibility with organic modifier.
- UV cut-off.
- Volatility for LC-MS.
- COLUMN STABILITY!

pH Range for Silica Columns



High pH Stability Comparison*




Exposure Conditions: Mobile phase, ACN/50mM Potassium phosphate buffer at indicated pH; Temperature, 30 °C.

LC Conditions: Mobile phase, ACN (or THF)/50mM Potassium phosphate buffer at indicated pH; Flow Rate, 1.0 mL/min.; Temperature, 30 °C; Injection Volume, 5 uL; Detection, 254nm.

* Column names are the trademarks of their respective manufacturers.

80 °C Aging Study

pH 7 ZirChrom-PBD vs. Xterra RP18*





LC Conditions: Mobile phase, ACN/50mM Potassium phosphate, pH 7.0; Flow rate, 1.0 ml/min.; Temperature, 80 °C, Injection Volume, 5 uL; Detection, 254nm. **25/75 ACN/Buffer ***30/70 ACN/Buffer



* Column names are the trademarks of their respective manufacturers.

ZirChrom[®]-PBD vs. Xterra[®] RP18

ZirChrom[®]-PBD







LC Conditions:
ZirChrom[®]-PBD; Mobile Phase, 28/72 acetonitrile/
20 mM potassium phosphate at pH=12.0; Flow Rate,
1.0 mL/min.; Temperature, 30°C; Detection, 254 nm.
Solute: 1=Labetalol, 2=Atenolol, 3=Acebutolol,
4=Metoprolol, 5=Oxprenolol, 6=Lidocaine,
7=Quinidine, 8=Alprenolol, 9=Propranolol.



LC Conditions:
Waters Xterra[™] RP18; Mobile Phase, 35/65 acetonitrile/
20 mM potassium phosphate at pH=12.0; Flow Rate,
1.0 mL/min.; Temperature, 30°C; Detection, 254 nm.
Solute: 1=Labetalol, 2=Atenolol, 3=Acebutolol,
4=Metoprolol, 5=Oxprenolol, 6=Lidocaine,
7=Quinidine, 8=Alprenolol, 9=Propranolol.

Xterra[®] RP18

ZirChrom[®]-PBD vs. Extend C18 Rapid Resolution



LC Conditions:
 ZirChrom[®]-PBD; Mobile Phase, 28/72 acetonitrile/
 20 mM potassium phosphate at pH=11.5; Flow Rate,
 1.0 mL/min.; Temperature, 40°C; Detection, 254 nm.
 Solutes: 1=Labetalol, 2=Atenolol, 3=Acebutolol,
 4=Metoprolol, 5=Oxprenolol, 6=Lidocaine,
 7=Quinidine, 8=Alprenolol, 9=Propranolol.

LC Conditions:
 Extend C18 Rapid; Mobile Phase, 45/55 acetonitrile/
 20 mM potassium phosphate at pH=11.5; Flow Rate,
 1.0 mL/min.; Temperature, 40°C; Detection, 254 nm.
 Solutes: 1=Labetalol, 2=Atenolol, 3=Acebutolol,
 4=Metoprolol, 5=Oxprenolol, 6=Lidocaine,
 7=Quinidine, 8=Alprenolol, 9=Propranolol.

The Preparation of Buffers with ZirChrom's Buffer Wizard

<http://www.zirchrom.com>

How to Prepare a Buffer

- Use Henderson-Hasselbalch equation to compute amount of acid and base needed for **simple monoprotic systems**.
- Get out your book on quantitative analysis and **calculate for an hour after studying for a day**.
- Use ZirChrom Buffer Wizard for all pH buffers.

<http://www.zirchrom.com>

Buffer Wizard - Microsoft Internet Explorer - [Working Offline]

File Edit View Favorites Tools Help

← Back → Home Search Favorites History

Address C:\Documents and Settings\administrator.BASSMASTER\Desktop\buffer_local.htm Go Links >>

ZirChrom Buffer Wizard™ Copyright© 2000 Aosheng Wang

Desired Buffer pH:	6.5	<ul style="list-style-type: none"> » Introduction to ZirChrom Buffer Wizard » How to Use ZirChrom Buffer Wizard » User-Defined Acid-Base » Key Issues in Buffer Preparation for HPLC » Do's and Don'ts of Using Buffers in HPLC » UV Cutoff Wavelength for Common pH Buffers » pH Ranges for Common HPLC Columns » Common Acid-Bases Used as pH Buffers in HPLC » ZirChrom Home
Acid:	<div style="border: 1px solid black; padding: 2px;"> Trifluoroacetic acid Acetic acid H3PO4 </div>	
Acid Stock Solution Conc. (M):	0.515 <input type="checkbox"/> 100% By wt%	
Desired Acid Conc. (M):	0.025	
Desired Total Buffer Volume (mL):	1500.0	
Base:	<div style="border: 1px solid black; padding: 2px;"> Tris Triethylamine Pyrrolidine </div>	
Base Stock Solution Conc. (M):	0.252 <input type="checkbox"/> 100% By wt%	
<input checked="" type="checkbox"/> Warning Calculate Report		
Required Acid Stock Volume (mL):	72.8155	
Required Base Stock Volume (mL):	173.5914	
% of Max Buffer Capacity:	55.5	

My Computer

Buffer Wizard Warning Panels

Microsoft Internet Explorer



Your buffer capacity may be too low at this pH. Suggest you choose a pH that is as close as possible to the pKa, choose a different buffer, or increase your buffer concentration. See Do's and Don'ts for additional information.

OK

Microsoft Internet Explorer



Only specially stabilized silica-based columns should be used at pH's $< 1.5-2$ or pH's $> 7-8$. The vast majority of all current RPLC phases are unstable under these conditions especially at temperatures above $40\text{ }^{\circ}\text{C}$. All silica-based RPLC phases are much less stable in phosphate and carbonate buffers at pH's > 7.5 and should not be used under these conditions. Zirconia based RPLC phases are stable at both very low and very high pH's, at high temperature ($150\text{ }^{\circ}\text{C}$) and in any buffer.

OK

Microsoft Internet Explorer



Your buffer concentration may be too high. Precipitation may result when the buffer is mixed with organic mobile phase modifier. See Do's and Don'ts for additional suggestions.

OK

Key Issues in Buffer Preparation for HPLC

- Adequate buffer capacity.
- Buffer solubility & compatibility with organic modifier.
- UV cut-off.
- Volatility for LC-MS.

The Big 14 Do's & Don'ts in Using Buffers in HPLC

1. Check pH of buffers *before* adding organic modifier not after.
2. Phosphate at $\text{pH} > 7$ destabilizes silica based columns. Should not be used > 40 °C. Silica columns are seriously destabilized at $\text{pH} < 2.5$ and $\text{pH} > 7-8$. *Use zirconia based columns for highest stability at any pH or temperature.*

Do's & Don'ts of Using Buffers in HPLC

3. Always *pre-check buffer solubility* by mixing the buffer with highest amount of organic modifier. Let sit for a minimum of 30 minutes to see if buffer precipitates. The buffer cation seriously effects its solubility. Sodium salts tend to be least soluble.

Do's & Don'ts of Using Buffers in HPLC

4. The buffer cation ($\text{Na}^+ < \text{K}^+ < \text{NH}_4^+ < \text{TEA}^+$) can help suppress silanophilic interactions in the order listed.
5. It is important to flush silica based columns with 15-20 volumes of buffer free eluent and then with pure methanol or pure acetonitrile (preferred) prior to overnight storage.
6. 200 ppm of sodium azide or alternatively a minimum of 20% (by volume) organic modifier should be added to the buffer to suppress bacterial growth.

Do's & Don'ts of Using Buffers in HPLC

7. 10-50 mM buffer is usually adequate. *25 mM is good to start* . Most inorganic buffers are more soluble in methanol than acetonitrile or THF.

Do's & Don'ts of Using Buffers in HPLC

8. If peak shape of acidic or basic analytes is poor increase the buffer capacity, or decrease the amount of sample injected. *Use more than 5 mM buffer and make sure that the buffer capacity is at least 20% of the maximum possible for the concentration used.*

Do's & Don'ts of Using Buffers in HPLC

9. Buffers containing ammonium or triethylammonium are more volatile than those containing Na^+ , K^+ , etc.
10. Buffers containing acetate, formate, fluoride, are more volatile than phosphate, citrate, etc.
11. Carbon dioxide can be lost by excessive sparging of the buffer. pH increases.

Do's & Don'ts of Using Buffers in HPLC

12. TFA and TEA degrade with time and the UV cutoff deteriorates.
13. TFA in glass ampoules gives better UV cutoff than bottled TFA.
14. Phosphate, has better UV cutoff than TFA, acetic acid or citric acid.

Effect of Organic Solvent on pH Buffers


- I. Canals, J. A. Portal, E. Bosch, M. Roses, “Retention of Ionizable Compounds on HPLC. 4. Mobile Phase pH Measurement in Methanol/Water”, *Anal. Chem.* 2000, 72, 1802-1809.

Thanks *very much*
for listening!


ZirChrom Separations & Cabot Corporation
Partners in Chromatography

PLEASE VISIT BOOTH 220-222

For more information and web access to the
free **Buffer Wizard**: **www.zirchrom.com**



ZirChrom Separations



SOLVENT	UV Cutoff (nm)	
	A=0.02	A=1.0
Methanol	250	205
Acetonitrile	210	190
THF	300	212
Isopropanol	250	205
Water	<190	<190

L(cm)	Column Dead Volumes (cc) ^a		
	Diameter (mm)		
	1	2.1	4.6
5	0.0271	0.119	0.573
10	0.0542	0.239	1.146
15	0.0812	0.358	1.719

Acid/Base	pKa
COOH	4.5-5.5
Aliph. N	9.0-11.0
Pyridine	5.17
Aniline	4.6

a. Total porosity = 0.69.

BUFFERS	pKa	Buffer Range	UV Cutoff (10mM) in nm
Trifluoroacetic*	0.5	1.5-2.5	210
H ₃ PO ₄	2.1	< 3.1	<200
	7.2	6.2-8.2	
	12.3	11.3-13.3	
Citric acid	3.1	2.1-6.4	230
	4.7		
	5.4		
Formic acid*	3.8	2.8-4.8	210
Acetic acid*	4.8	3.8-5.8	210
Carbonate	6.4	5.4-7.4	<210
	10.3	9.3-10.3	
Bis-tris propane	6.8	5.8-7.8	215
	9.0	8.0-10.0	
Tris	8.3	7.3-9.3	205
Ammonia*	9.2	8.2-10.2	200
1-methylpiperidine	10.1	9.1-11.1	215
Triethylamine*	11.0	10.0-12.0	<200

* denotes volatile buffer

Solvent	Acetic acid	Acetone	Acetonitrile	Carbon tetrachloride	Chloroform	Cyclohexane	Cyclopentane	Dichloromethane	Dimethylformamide	Dioxane	Ethyl acetate	Ethanol	Diethyl ether	Hexane	Methanol	Pentane	n-Propanol	Diisopropyl ether	Tetrahydrofuran	Toluene	Water	
Acetic acid																						
Acetone																						
Acetonitrile																						
Carbon tetrachloride																						
Chloroform																						
Cyclohexane																						
Cyclopentane																						
Dichloromethane																						
Dimethylformamide																						
Dioxane																						
Ethyl acetate																						
Ethanol																						
Diethyl ether																						
Hexane																						
Methanol																						
Pentane																						
n-Propanol																						
Diisopropyl ether																						
Tetrahydrofuran																						
Toluene																						
Water																						

Legend: ■ Immiscible, Miscible

To access our FREE automated HPLC Buffer Wizard calculator and a great deal of other HPLC data visit our web site at <http://www.zirchrom.com>

World's Most Robust HPLC Columns. Telephone: 763-421-5264, Fax: 763-421-2319



Product List

Part #	Product Name	Chromatographic Mode and Uses
ZR01	ZirChrom [®] -CARB	Reversed-Phase
ZR02	ZirChrom [®] -PHASE	Normal Phase and SEC
ZR03	ZirChrom [®] -PBD	Reversed-Phase
ZR04	ZirChrom [®] -WCX	Weak Cation-Exchanger
ZR05	ZirChrom [®] -WAX	Weak Anion-Exchanger and Sugar Analysis
ZR06	ZirChrom [®] -SAX	Strong Anion-Exchanger
ZR07	ZirChrom [®] -SHAX	Strong Hydrophilic Anion-Exchanger
ZR08	ZirChrom [®] -PEZ	Cation-Exchanger for Proteins