



Zirconia: the Ideal Substrate for Ion-Exchange LC and LC-MS

EAS 2005

Bingwen Yan¹, Clayton V. McNeff¹, Richard A. Henry², and David S. Bell²

¹ZirChrom Separations, Inc., 617 Pierce Street, Anoka, MN 55303

²Supelco, 595 North Harrison Road, Bellefonte, PA 16823

ZirChrom

1-866-STABLE-1
www.zirchrom.com

... For Peak Performance




Ion-Exchange LC (HPLC)

- **Most ion-exchange LC has been done with polymer-based support particles.**
 - Polymers often exhibit poor efficiency and are subject to dimensional change when exposed to high temperature, high pressure, high flow rate and organic solvents.
- **Silica has not been highly successful as a base for ion-exchange.**
 - Substrate has weak cation exchange properties which are not easy to reproduce and control; silanols overload easily and reduce column performance.
 - Coated and bonded phases can exhibit limited stability in aqueous solution over range of pH and temperature conditions needed for ion-exchange.
- **Zirconia may be an ideal particle for ion-exchange.**
 - Zirconia is rigid and very stable over a wide range of operating conditions.
 - Both cation and anion exchange versions are readily prepared.
 - Rugged industrial and biochemical ion separations may be developed or transferred by following standard rules for IEC or IC method development.
- **Challenge- can ion-exchange LC-MS become a routine tool?**



Review: Zirconia Chemistry




- Ligand exchange interaction : $\text{Zr-L} + \text{A}^- \ll \text{Zr-A}^- + \text{L}$
- When ligand is charged, surface becomes charged!

* Base is electron donor; acid is electron acceptor; more general than Brönsted.




Four Zirconia-based Options for IEC



- **Bare zirconia**
- Phosphate, fluoride and other anionic additives that are replaceable.
- primarily a strong cation exchanger (SCX)



- **Zirconia with PEI coating**
- Cross-linked to resist removal even under extreme conditions
- Weak anion exchanger (WAX)



- **Zirconia with EDTPA chelator modification**
- Multidentate attachment
- Very stable, but can be replaced or restored
- Strong cation exchanger (SCX)




- **Zirconia with quaternized PEI coating**
- Cross-linked to resist removal even under extreme conditions
- Strong anion exchanger (SAX)

Time (min)



Proteins on Zr-PO₄ and ZirChrom-PEZ (EDTPA)




- 1. Myoglobin
- 2. Ribonuclease
- 3. Chymotrypsin
- 4. Lysozyme
- 5. Cytochrome *c*

Mobile phase, 50 to 500 mM K₂HPO₄ in 20 min., pH 7.0; Flow Rate, 1 mL/min; Detection, 280 nm.



Proteins on ZirChrom-PEZ




Gradient elution of cytochrome *c* from horse heart. Mobile phase, 4 mM EDTPA, 20 mM MES [2-(N-morpholino)ethanesulfonic acid], pH 5.5 with a linear gradient of 0 to 1 M NaCl in 30 min; Detection, 410 nm.



ZirChrom®

Small-Scale Purification of MAB IgG_{2a} Contaminated with BSA




Note: BSA is not retained under sample loading conditions.

LC Conditions: ZirChrom-PEZ, 5 x 4.6mm (part#ZR08-0546); 100 μ l injection of BSA (6.0 mg/ml) contaminated MAB (1.0 mg/ml) eluted by salt step gradient; Mobile phase, 20 mM MES, 4 mM EDTPA, 0.05 M-to-1.0 M NaCl pH=5.5; Flow rate, 2.0 ml/min; Temperature: 30°C; Detection, 280 nm.



Zirconia-PEI (Zr-PEI) for Anion Exchange


- Coat and crosslink polyethyleneimine (PEI) for weak anion exchange (WAX).
- Quaternize with methyl iodide for strong anion exchange (SAX).



Lewis acid sites can compete for anions unless blocked by a stronger Lewis base.



Small Anions on ZirChrom -WAX




Addition of organic solvent may improve peak shape of organic anions by eliminating any RP effect.

LC Conditions:

Column, ZirChrom®-WAX, 150 x 4.6 mm i.d. (part# ZR05-1546); Mobile Phase, 45mM ammonium phosphate dibasic at pH 8.2; Flow Rate, 1.0 ml/min; Detection, 240 nm; Column Temperature = 40°C. Solutes: 1 = bromate, 2 = nitrite, 3 = benzoic acid, 4 = nitrate, 5 = p-chlorobenzoic acid, 6 = p-bromobenzoic acid, 7 = iodide, 8 = p-fluorobenzoic acid, 9 = p-iodobenzoic acid.



Trace Iodide Separation on ZirChrom®-SAX




Analytes:
2M nitrate
sample matrix
and Iodide

LC Conditions:

Column, ZirChrom®-SAX, 50 x 4.6 mm, (part# ZR06-0546); Mobile Phase: 25mM ammonium phosphate, 275mM NaCl at pH 8.0; Flow Rate, 1.0 mL/min; Temperature, 30°C; Detection, 226 nm; Injection volume: 5 ul.



“Green” Analysis of Diet Soft Drinks Containing Caffeine and Aspartame




LC Conditions:

Column, ZirChrom®-SAX, 100 x 3.0 mm (part# ZR06-1030); Mobile Phase, 10mM Ammonium phosphate, 5mM Ammonium carbonate, pH 6.6; Flow rate, 1.0 ml/min; Temperature, 50 °C; Injection Vol., 5.0 ml; Pressure Drop, 205 bar; Detection, 210 nm.



LC/MS Compatible Separation of Non-Steroidal Anti-Inflammatory Drugs (NSAIDs) on ZirChrom[®]-SAX



1-Acetaminophen

2-Ibuprofen

3-Naproxen

4-Impurity


5-Ketoprofen

LC Conditions:

Column, ZirChrom[®]-SAX, 50 x 4.6 mm i.d. (part number: ZR06-0546); Mobile Phase, 80/20 ACN/15 mM ammonium formate, pH=4.0 (adjusted with formic acid); Flow Rate, 1.0 ml/min; Temperature, 35 °C; Injection Vol., 1.0 ml; Detection, 254 nm.



Water-Soluble Vitamin Analysis on ZirChrom[®]-SAX



- 1 - Thiamine (Vit. B₁)
- 2 - Pyridoxine (Vit. B₆)
- 3 - Nicotinamide (form of Vit. B₃)
- 4 - Riboflavin (Vit. B₂)
- 5 - Nicotinic acid (form of Vit. B₃),
- 6 - Ascorbic acid (Vit. C)

**Vitamin C is largely retained
on ZirChrom-SAX**

LC Conditions:

Column: ZirChrom[®]-SAX, 150 x 4.6 mm i.d. (part number: ZR06-1546); Mobile Phase, 50 mM Ammonium dihydrogenphosphate, pH 4.5; Flow rate, 1.0 ml/min; Temperature, 30 °C; Injection Vol., 5.0 ml; Detection, 254 nm.



ZirChrom®

Conclusions

- Zirconia is a very attractive LC substrate because of its unparalleled mechanical and chemical stability (especially in aqueous solution).
- It has great potential to become a primary substrate for both preparative and analytical ion-exchange LC with UV, conductivity, electrochemical and other common detectors that can tolerate phosphate, fluoride, chloride, sulfate and other nonvolatile additives.
 - Retention rules for ions are easily understood.
 - Ionic strength gradients are tolerated by the packing (chlorides, sulfates, nitrates, etc.).
 - Zirconia and the IEC mechanism are both tolerant of organic solvent to counter excessive hydrophobic solute retention or to elute hydrophobic interferences at the solvent front.
- Zirconia is potentially very useful for ion-exchange LC-MS of both positive ions (SCX mode) and negative ions (WAX and SAX modes) using volatile mobile phases or online clean-up.
 - Ammonium acetate, ammonium formate, ammonium carbonate and ammonium hydroxide additives should be useful; effect of additives and ionic strength gradients on MS signal requires further investigation.
 - IEC mode tolerates high organic to reduce RP effects and maximize MS-ESI signal.



For More Information, Visit Us at Booth 220

A white rectangular banner containing the ZirChrom logo on the left, the text 'ZirChrom®' in the center, and a blue chromatogram on the right. The chromatogram shows a baseline with three distinct peaks of varying heights. To the right of the peaks, the text '1-866-STABLE-1' and 'www.zirchrom.com' is displayed. Below the chromatogram, the tagline '... For Peak Performance' is written in a smaller blue font.

ZirChrom®

1-866-STABLE-1
www.zirchrom.com

... For Peak Performance

For more information and web access to the free
Buffer Wizard: www.zirchrom.com