



Evaluation of 1.7 μ m Porous Zirconia Particles for UHPLC

Bingwen Yan¹, Clayton V. McNeff¹, R.A. Henry², Dan Nowlan¹

1. ZirChrom Separations, Inc. 617 Pierce St., Anoka, MN 55303,

2. Independent Consultant, 983 Greenbriar Drive, State College, PA 16801




Outline

- **Chemistry of zirconia-based HPLC columns**
 - **Surface chemistry**
 - **Eluotropic series**
 - **Polymer coated zirconia**
- **Zirconia phases for Multi-Mode applications and unique selectivity**
 - **Zr-PBD and Zr-PS [RP and SCX in phosphate mobile phases]**
 - **Zr-MS [RP and SCX in any mobile phase (MS friendly)]**
 - **Zr SAX [RP and SAX in any mobile phase]**
- **Evaluation of sub-2 μ m Zr-PBD for UHPLC**
 - **Specifications**
 - **Performance (van Deemter Plots)**
 - **Applications**



Surface Chemistry of Zirconia




Zirconia chemistry is dominated by Lewis acid-base reactions (electron donor-acceptor) involving empty d-orbitals.



Other Lewis base examples: PO_4^{3-} , RCO_2^- , F^- , Catechol, etc.




Interaction Strength of Lewis Bases with Zirconia¹

Interaction Strength	Lewis Base (L)
Strongest	Hydroxide
	Phosphate
	Fluoride
	Citrate
	Sulfate
	Acetate
	Formate
	Nitrate
	Chloride
	Water
	Weakest

Small Lewis base (anions) with high electron density and low polarizability interact more strongly with Zr atoms.

Strong Lewis base solutes adsorb to Zr atoms unless displaced by stronger bases in mobile phase.




Multi-Mode Behavior of Zirconia

- Zirconia substrate exhibits reproducible polar and ionic solute interaction.
- With stable organic coatings, reproducible reversed-phase behavior can be added.
- Extreme resistance to temperature, pH and mechanical stress are unique advantages.



Addition of RP Behavior in Coated Zirconia Phases²



- Retention (and selectivity) of ionic analytes modulated by pH, buffer/salt type and concentrations, and temperature.
- Retention of neutral solutes modulated by organic solvent.



ZirChrom®


Definition of Multi-Mode HPLC

- Certain columns have two or more HPLC retention modes working simultaneously.
- Can be called Multi-Mode, Dual-Mode or Mixed-Mode.
- Typically, modes include **Reversed-Phase**, **Ion-Exchange** (anion or cation) and **Normal-Phase** (or **HILIC** in aqueous mobile phases).
- These modes are orthogonal to one another (analytes usually elute in different order).



C18 and Zr-PBD are Orthogonal for Basic Drugs²

C18 (RP) columns separate mainly by hydrophobic forces and Zr-PBD columns separate by a combination of ionic and hydrophobic forces




Zr-PBD and Si-C18 have very different selectivity for ionic drugs under phosphate conditions due to the SCX component of ZrO_2 .

Solutes

1. Chlordiazepoxide
2. Hydroxyzine
3. Buclizine
4. Thiothixene
5. Doxepin
6. Amitriptyline
7. Imipramine
8. Perphenazine
9. Nortriptyline
10. Desipramine
11. Thioridazine



IEX Mode Increases Retention for Zr-PBD vs Si-C18²




LC Conditions: Machine-mixed 80/20 ACN/10 mM ammonium acetate pH=6.7 without pH adjustment; Flow rate, 1.0 mL/min.; Injection volume 0.1 μ L; Temperature, 35 $^{\circ}$ C; Detection at 254 nm; Columns, Zr-PBD, 50 x 4.6 mm i.d. (3 μ m particles); Silica-C18 150 x 4.6 mm i.d., (3.5 μ m particles).



Ion-Exchange Mode on Zirconia³


Quaternary amines paraquat and diquat are retained and resolved on Zr-PS (also Zr-PBD or bare ZrO₂) due to a **cation exchange** mechanism; 50% ACN is useful to suppress phase retention by RP mode.

Silica-C18:
reversed-phase (silanol tailing)



column: Discovery C18, 15 cm x 4.6 mm I.D., 3 μ m
mobile phase: 5% acetonitrile in 25 mM phosphate (pH 7)
flow rate: 1 mL/min.
temp.: 35 °C
det.: UV 290 nm

Zirconia-PS:
primarily ion-exchange




column: Discovery Zr-PS, 7.5 cm x 4.6 mm, 3 μ m
mobile phase: 50% acetonitrile in 25 mM phosphate (pH 7)
flow rate: 3 mL/min.
temp.: 65 °C
det.: UV 290 nm



Anticholinergics on Zr-PBD³

Quaternary amines (except 2)



LC Conditions

Discovery Zr-PBD, 100mm x 2.1mm i.d., 3 μ m


Mobile Phase A: 50:50 [20 mM H₃PO₄, pH 7.0 w/ NH₄OH]:water

Mobile Phase B: 50:30:20 [20 mM H₃PO₄, pH 7.0 w/ NH₄OH]:water:ACN


Gradient: 10-100% B over 18 minutes

Temp: 80 °C, Flow: 0.3 mL/min, Inj vol: 2 μ L in 60% MeOH,


UV@225 nm




1, Pipenzolate (20 mg/L)




2, Scopolamine (100 mg/L)




3, Ipratropium (100 mg/L)



4, Methscopolamine (100 mg/L)



5, Propantheline (20 mg/L)




6, Oxyphenonium (100 mg/L)




Surface Chemistry and Retention Mechanisms of QPEI-Zirconia⁴

- Anion-exchange
- Hydrophobic interactions
- Lewis acid-base interactions





Water-Soluble Vitamin Analysis on ZirChrom[®]-SAX⁴



- 1 - Thiamine (Vit. B₁)
- 2 - Pyridoxine (Vit. B₆)
- 3 - Nicotinamide (form of Vit. B₃)
- 4 - Riboflavin (Vit. B₂)
- 5 - Nicotinic acid (form of Vit. B₃)
- 6 - Ascorbic acid (Vit. C)

→ Acid anions well retained on SAX

LC Conditions: Column: ZirChrom[®]-SAX, 150 x 4.6 mm i.d. (part number: ZR06-1546),
Mobile Phase: 50 mM Ammonium dihydrogenphosphate, pH 4.5, Flow rate: 1.0 mL/min.
Temperature: 30 °C, Injection Vol.: 5.0 µL, Detection: UV at 254 nm





Evaluation of sub-2 μ m Porous Zirconia Particles for UHPLC

- Since about 1970, there has been a steady effort to develop smaller particle particles that improve HPLC column efficiency. Most recently, columns with sub-2 μ m particles have been given the name Ultra-HPLC (UHPLC).
- The performance of a sub-2 μ m zirconia particle coated with PBD has been investigated.



Analytical Diameter Porous Zirconia Particles

Particles are produced by a sol-gel process using ca. 1000Å sol.



1 μm 25000X




Properties of Porous Analytical Zirconia

<u>Characteristic</u>	<u>Property</u>
Surface area (m ² /g)	22
Pore volume (cc/g)	0.13
Pore diameter (Å)	250-300
Porosity	0.45
Density (gm/cc)	2.6 (2.5x silica)
Particle diameter (μm)	3.0 (and sub-2)



ZirChrom®

Sub-2 μm Particle Size Distribution Compared to 3 μm



Volume Statistics (Arithmetic)


9573-5C.01

Volume:	100%		
Mean:	1.923 μm	95% Conf. Limits:	1.919-1.927 μm
Median:	1.919 μm	S.D.:	0.426 μm
Mean/Median ratio:	1.002	Variance:	0.182 μm^2
Mode:	1.984 μm	C.V.:	22.2%



Sub-2 μ m BET Data*

Adsorption Pore Volume Plot, ($dV/d\log D$)




*Data identical to 3 μ m zirconia

Pittcon 2009



Idealized van Deemter Plots


$$H = A + B/v + Cv$$



Elements of the drawing provided by Sigma-Supelco



Comparison Between Sub- $2\mu\text{m}$ and $3\mu\text{m}$ Particles: Alkylbenzenes



Columns: ZirChrom PBD, 50 x 4.6mm; Mobile phase: 50/50 ACN/water; Flow 2.0 mL/min (0.53 cm/sec); Temperature: 30 °C; UV@254nm; Agilent 1100 with micro flow cell.



Flow Studies on 3 μ m Zr-PBD: Alkylbenzenes




Plate height based on van Deemter Equation vs linear velocity at various temperatures for retained solutes: Alkylbenzenes, Temperature: 30 °C, Mobile phase: 55/45 ACN/water, Column: ZirChrom PBD, 50 x 4.6mm, Agilent 1100/UV with Micro Cell (0.007" i.d. tubing).



Flow Studies on Sub-2 μm Zr-PBD: Alkylbenzenes





Plate height based on van Deemter Equation vs linear velocity for retained solutes: Alkylbenzenes, Temperature 30 °C, Mobile phase: **50/50 ACN/water** (keep k' in the same range as 3 μm particles), Column: 50 x 4.6mm, Agilent 1100/UV with Micro Cell (0.007'' i.d. tubing).





Comparison Between Sub-2 μm and 3 μm Zr-PBD



Resistance to mass transfer is **reduced** as the particle size is decreased. Fast analysis can be achieved at higher flow rate with minimal efficiency loss.



Zr-PBD Separation of Catecholamines: **Sub-2 μ m** vs **3 μ m**



Conditions

Column: ZirChrom-PBD 50 x 4.6mm
Mobile Phase: 85/15 ACN/30mM NH₄OAc,
10mM NH₄H₂PO₄ adjusted to
pH=3.4 w/ HCl

Flow rate = 1.5 mL/min

Temperature = 30°C


Inj Vol = 5 μ L

Elution order: 1=Tyramine, 2=Epinephrine
3=Dopamine,
4=3,4-dihydroxynorephedrine

- **Smaller particle shows increased efficiency**
- **Identical selectivity to 3 μ m**
- **Separation time can readily be decreased to <1min for sub-2 μ m columns.**





Sub-2 μ m Pressure Drop at Different Temperatures





UHPLC Column Hardware



- Column hardware must withstand very high packing and operating pressures.
- Frits must contain particles in the $0.5\mu\text{m}$ range or smaller.
- Composite frits are composed of $0.5\mu\text{m}$ glass fiber mats sandwiched between $2\text{-}3\mu\text{m}$ stainless wire meshes.



Summary and Conclusions

- Multi-mode HPLC columns have become popular for difficult applications where compounds have ionic character and vary widely in chemical nature.
- Several ZirChrom phases are ideal and popular for multi-mode applications, including Zr-PBD, Zr-PS, Zr-MS and Zr-SAX, and are stable over a much wider range of pH and temperature than any silica-based phase.
- ZirChrom phases are currently available in high efficiency 3 μ m analytical particles.
- A new sub-2 μ m particle has been developed for preparing UHPLC columns that produce very high efficiency in excess of 200,000 N/M with PBD polymer-coated phase. Additional phase coatings are being investigated.



References and Acknowledgements

1. **J. A. Blackwell and P. W. Carr, "Development of an Eluotropic Series for the Chromatography of Lewis Bases on Zirconium Oxide," Anal. Chem. 64, 863-73 (1992).**
2. **R. A. Henry, H. K. Brandes, D. S. Bell and C. T. Santasania, 30th Annual HPLC Meeting, Oral Presentation, 2006, San Francisco, CA.**
3. **R. A. Henry and H. K. Brandes, Eastern Analytical Symposium, Oral Presentation, 2006, Somerset, NJ.**
4. **B. Yan. C. V. McNeff. R. A. Henry and D. Nowlan, Eastern Analytical Symposium, Poster, 2008, Somerset, NJ.**

The assistance of Hillel Brandes of Supelco and Kelly Johnson of ZirChrom in providing data and other helpful comments is greatly appreciated.