



High Temperature Liquid Chromatography

C.V. McNeff^{1*}, B. Yan¹, D. R. Stoll², R.A. Henry³

¹ ZirChrom Separations, Inc., 617 Pierce Street, Anoka, MN 55303

**² Department of Chemistry, University of Minnesota, 207 Pleasant St.
Minneapolis, MN 55455.**

³ Independent Consultant, 983 Greenbriar Drive, State College, PA 16801



OUTLINE

- **Advantages of High Temperature HPLC**
 - Theoretical Effects of High Temperature HPLC
 - Practical Analytical Advantages of Using High Temperature HPLC
- **Development of New Stationary Phases**
 - Selectivity Comparison of Zirconia Based Stationary Phases with C18 Silica and Other Columns
 - High Temperature Separations
- **Using Temperature to Control Selectivity**
 - Importance of Selectivity in HPLC Optimization




Theoretical Advantages to High Temperature LC

van Deemter Plot

$$h = A + \frac{B}{v} + C v + D v^{2/3} + \frac{3D_m}{8k_d d_p^2} v$$

R. D. Antia and Cs. Horvath, *J. Chromatogr.*, 435, 1-15 (1988).

**Practical Limit
Temperature Dependence**

$$\frac{t}{N} \propto (1 + k') \frac{L^{2/3}}{\Delta P_{\max}^{2/3}} \frac{\eta}{T^{1/3}}$$



Guiochon, Georges, *Anal. Chem.*, 52, 2002-2008 (1980).

Three ways that temperature increases efficiency and speed

- Increased temperature increases diffusivity, thus decreasing the reduced velocity
- Increased temperature accelerates sorption kinetics
- Increased temperature decreases mobile phase viscosity




Theoretical Effect of Temperature on Column Efficiency






Estimated Effect of Temperature on Viscosity*



*H. Chen and Cs. Horvath, "Rapid Separation of Proteins by RP-HPLC at Elevated Temperatures," *Anal. Methods Instrum.*, **1**, 213-222 (1993).




Effect of Temperature on Theoretical Analysis Time at Constant Pressure and Plate Count*



*R. D. Antia and Cs. Horvath, *J. Chromatogr.*, **435**, 1-15 (1988).



Practical Advantages of Column Stability






List of HTLC compatible reversed-phase columns

Manufacturer	Column Name	Stationary Phase Type	Temperature Limit (°C)	Selectivity versus C18
Polymer Laboratories	PLRP	Polymer	200	Different
Selerity	Blaze	Silica	200	Different
Supelco	DiscoveryZR-Carbon	Carbon Clad Zirconia	200	Different
Supelco	DiscoveryZR-CarbonC18	Modified Carbon on Zirconia	200	Different
Thermo-Electron	Hypercarb	Carbon	200	Similar
ZirChrom Separations, Inc.	ZirChrom-CARB	Carbon Clad Zirconia	200	Different
ZirChrom Separations, Inc.	Diamondbond-C18	Modified Carbon on Zirconia	200	Different
Jordi	Jordi DVB	Polymer	150	Different
Sachtleben	Sachtopore-RP	Polymer coated Titania	150	Different
Supelco	DiscoveryZR-PBD	Polymer Coated Zirconia	150	Different
Supelco	DiscoveryZR-PS	Polymer Coated Zirconia	150	Different
ZirChrom Separations, Inc.	ZirChrom-PS	Polymer Coated Zirconia	150	Different
ZirChrom Separations, Inc.	ZirChrom-PBD	Polymer Coated Zirconia	150	Different
Agilent	SB Extend-C18	Silica	90	Similar
Waters	X-Bridge	Silica	80	Similar



Stationary Phase Comparison

Average Scatter of κ - κ Plots for Two Kinds of Stationary Phases Using 22 Solutes




- Carbon-ZrO₂
- PBD-ZrO₂
- C18-SiO₂ (ODS)
- Phenyl-SiO₂
- CN-SiO₂
- PRP

➤ For *non-electrolytes*, C-ZrO₂ and aliphatic phases have the most different selectivities.



Fast Separations Non-Steroidal Anti-Inflammatories




LC Conditions: Column, 50 x 4.6 DiamondBond™-C18; Mobile phase, 25/75 ACN/40mM phosphoric acid, pH 2.3; Flow rate, 5.5 ml/min.; Temperature, 150 °C; Injection volume, 1ul; Detection at 254nm; Solute concentration, 0.15 mg/ml.; Solutes, 1= Acetaminophen, 2=Ketoprofen, 3=Naproxen, 4=Ibuprofen, 5=Oxaprofen.



ZirChrom®

Fast β -Blockers Separation

Column Temperature = 150°C, pH = 11



LC Conditions: Column, 50 x 4.6 Diamondbond-C18, OD0121601A; Mobile phase, 45/55 ACN/20mM Ammonium Phosphate pH11.0; Flow rate, 3.0 ml/min; Temperature, 150 °C; Injection volume, 1.0 ul; Detection at 210 nm; Solutes, 1=Labetalol, 2=Metoprolol, 3=Alprenolol




Resolution: The Importance of Selectivity

Efficiency	Retention	Selectivity
↓	↓	↓
$R = \frac{\sqrt{N}}{4}$	$\frac{k'}{k'+1}$	$\frac{\alpha-1}{\alpha}$

$$\alpha = \frac{k_j'}{k_i'}$$

➤ Selectivity (α) has the greatest impact on improving resolution.







ZirChrom®

Comparison of Variables Affecting Selectivity

30% ACN vs. 50% ACN




80°C vs. 30°C




Stationary Phase Type


Carbon-ZrO₂ vs.
PBD-ZrO₂



MeOH vs. THF



C18-SiO₂ vs. PBD-ZrO₂




- Stationary phase type has a very large effect on selectivity.



ZirChrom®

Thermally Tuned Tandem Columns (T³C)


A Mechanism to Continuously Adjust the Stationary Phase




Column 1



Column 2




**Optimized
T³C**





Separation of Ten Triazine Herbicides by T³C




Solutes:

- | | |
|--------------|------------------|
| 1. Simazine | 6. Ametryn |
| 2. Cyanazine | 7. Propazine |
| 3. Simetryn | 8. Terbutylazine |
| 4. Atrazine | 9. Prometryn |
| 5. Prometon | 10. Terbutryn |

Other conditions:

30/70 ACN/water

1ml/min; 254 nm detection



➤ T³C can improve separation without increasing analysis time.



Conclusions

- (1) Zirconia Based Stationary Phases are *ultra-durable* and *efficient*, *stable* at the *extremes of pH* and at column temperatures as high as *200°C*.
- (2) *ZirChrom®-CARB* has the *most different selectivity* relative to conventional ODS phases for the 22 selected non-ionizable compounds.
- (3) High Temperature Liquid Chromatography (HTLC) is a *powerful technique* that can be used as a *routine analytical tool* in the development of separation methods.
- (4) HTLC is a *unique tool* in altering chromatographic selectivity (*T³C method*), increasing analysis speed.
- (5) HTLC capability will become an *important* part of HPLC *system design* in order to fully utilize the benefits of columns prepared with ultra-small particles and ultrafast analyses.