


Polymer Coated Titanium Dioxide as a Reverse Phase Sorbent in Liquid Chromatography



**Presented at the PREP 2004 Conference
May 23-26 in Baltimore/USA**

**J. Winkler ^{a*}, B. Hirthe ^a, G. Benkner ^a, J. Blome ^b,
C. Spitzwieser ^a, R. Wittenberg ^a**

^aSachtleben Chemie GmbH, Duisburg, Germany

^bHPLC-Analytik Johannes Blome, Erolzheim

*to whom correspondence should be addressed

Contents

Chapter	Topic	Page
1	Titel	1
2	Contents	2
3	Introduction	3 - 4
4	List of Separations	5
5	Alphabetical Listing of drugs	6
6	RP - Separations	7 - 14



■ Introduction

Porous titanium dioxide is different from other inorganic sorbents used in high pressure liquid chromatography in the sense that it has a high resistance towards the attack of acids as well as bases. Titanium dioxide in the anatase or rutile crystal modification is furthermore completely inert against moisture uptake. For that reason polymer modified titanium dioxide is useful as a sorbent for HPLC under extreme conditions. RP titanium dioxide from Sachtleben carries the name Sachtopore RP .

■ Physical Properties of Sachtopore

Sachtopore	100 Å	200 Å	330 Å	2000 Å
Pore Diameter	100	200	300	2000
Specific surface area [m ² /g] (typical values)	55	30	15	<5
Specific pore volume [ml/g] (typical values)	0,21	0,17	0,12	n.d.
Mechanical stability				
Applicable packing pressure [bar] (typical values)	750	750	750	600

■ Particle Sizes

Sachtopore is available in particle sizes of 3 µm, 5 µm, 10 µm, 20 µm, 40 µm and 80 µm.

■ Samples and Elution Technique

The samples were chosen under the aspect of a high or low pKa value, respectively, and the lack of acidic next to basic groups in the molecule which would lead to ionic structures due to intramolecular acid-base reactions. To suppress dissociation basic samples were run under alkaline conditions, whereas acidic samples were run at low mobile phase pH values. All separations were carried out in an isocratic mode.

■ pH-range

Mobile phases with pH values between 10 and 1,5 were used. Borate buffers were tested because of their high pH values and their low UV cut-off wavelength of 210 nm, making them interesting when UV detectors are used.

■ Results

From the chromatograms the plate numbers per mm and peak asymmetries for each substance were determined. The quotient of plate number to asymmetry was taken as a measure for column performance:

	pH	Particle size [µm]	Average Plate Number [mm-1]	Average asymmetry	Performance plate number asymmetry
Fig. 1	10	3	43,1	1,15	37,5
Fig. 2	10	3	44,1	1,27	34,7
Fig. 3	10	5	28,7	1,54	18,6
Fig. 4	10	5	28,3	1,27	22,3
Fig. 5	9,5	5	21,0	1,06	19,8
Fig. 6	9,5	5	22,8	1,31	17,4
Fig. 7	9,5	5	19,8	1,14	17,7
Fig. 8	9	5	35,3	1,38	25,6
Fig. 9	9	3	42,6	1,27	33,5
Fig. 10	2,3	3	47,1	1,13	41,7
Fig. 11	2,3	3	45,2	1,31	34,5
Fig. 12	2,2	5	31,8	1,16	27,4
Fig. 14	1,5	3	45,2	1,30	34,8
Fig. 15	1,5	5	29,7	1,24	23,9
Fig. 16	1,5	5	27,9	1,27	22,0

■ Preparative LC simulation

To simulate preparative chromatography conditions, Pentifylline was extracted from a tablet of Cosaldon™. Using 10 mM Na₂B₄O₇ + 1 mM H₃BO₃ / 20% Acetonitrile (pH 9 in aqueous component) as mobile phase the chromatogram shown in Fig. 13 was obtained under overload conditions. Impurities were separated at retention time of 1.55, 2.07 and 4.22 minutes. The Pentifylline peak appeared with an almost rectangular shape.

■ Conclusion

Polymer modified Sachtapore is useful as a reverse phase sorbent in the separation for samples with high acidities or basicities and/or under extreme pH values of the mobile phase. The use of reverse phase Sachtapore in preparative chromatography under similar challenging conditions appears promising.

List of Separations


Type of Separation	Fig.-Number	Page
Separation of Antipsychotic Drugs	Fig. 1	1
Separation of Antipsychotic Drugs	Fig. 2	1
Separation of Antihistamines	Fig. 3	2
Separation of Cardiac Drugs	Fig. 4	2
Separation of Indole Alkaloids	Fig. 5	3
Separation of Antihypertensives	Fig. 6	3
Separation of Ergot Alkaloids	Fig. 7	4
Separation of Antipyretic Drugs	Fig. 8	4
Separation of Antihypertensives	Fig. 9	5
Separation of Antipsychotic Drugs	Fig. 10	5
Separation of Antipsychotic Drugs	Fig. 11	6
Separation of Antipsychotic Drugs	Fig. 12	6
Preparative Separation of Pentifyline	Fig. 13	7
Separation of Antiphlogistic Drugs	Fig. 14	7
Separation of Antiseptics	Fig. 15	8
Separation of Antiseptics	Fig. 16	8

Alphabetical Listing of Drugs

Drug	in Fig.	Drug	in Fig.
Aminoantipyrine	8	Isosorbide dinitrate	6 / 9
Aminophenazone	8	Olanzapine	2
Amiodarone	4	Oxyphenylbutazone	11/12/14
Antipyrine	8	Paracetamol	8
benzoic acid	15	Pentifylline	9/13
Betahistine	6	Phenazine	12
Carbinoxamine	3	Phenol	15
Chlomipramin	1	Phenylbutazone	11/12/14
Chlorphenamine	3	Phenyltoloxamine	3
Clofibrate	6	p-nitro benzoic acid	15
Clopazine	1	Paracetamole	8
Clozapine	2	Prazosin	6 / 9
Clomipramine	1	maleic acid	3
Diclofenac	10	Medazepame	2
Diflunisal	10/11/12	Metoprolol	4
Dihydroergocristine	7	Naproxen	11/12/14
Dihydroergocornine	7	Nifedipine	4
Dihydroergocryptine	7	Promethazine	1
Dihydroergotamin methylsulfate	6	Propylphenazone	11/12/13
Diphenylpyraline	3	Reserpine	5
Dipyridamole	4	Rofecoxib	12
Doxepin	2	salicylic acid	15
Ergotamine	6	Thymol	16
Etafenone	4	Trapidil	4
ethenzamide	8	Verapamil	4
Flupentixol	1	Vincamin	9
fumaric acid	15	4-hydroxybenzoic acid ethyl ester	16
Haloperidol	1	4-hydroxybenzoic acid methyl ester	16
Ibuprofene	11/14	4-hydroxybenzoic acid propyl ester	16
Indometacine	10/11		

Fig. 1: RP Separation of Antipsychotic Drugs

pH 10



Column

Sachtopore RP (300 Å, 4 µm, 150 :

Eluent

72% Water
(+ 20 mM Na₂CO₃ + 10 mM NaHC
28% Acetonitrile

Flow


1 ml/min

Detection

UV, 250 nm

Fig. 2: RP Separation of Antipsychotic Drugs

pH 10



Column

Sachtopore RP (300 Å, 3 µm, 150 x 4 mm)

Eluent

72% Water
(+ 20 mM Na₂CO₃ + 10 mM NaHCO₃)
28% Acetonitrile

Flow


1 ml/min

Detection

UV, 250 nm

Fig. 3: RP Separation of Antihistimines

pH 10



Column

Sachtopore RP (300 Å, 3 µm, 150 x 4 mm)

Eluent

70% Water
(+ 10 mM Borax + 10 mM Soda)
30% Acetonitrile

Flow


1 ml/min

Detection

UV, 230 nm

Fig. 4: RP Separation of Cardiac Drugs

pH 10



Column

Sachtopore RP (300 Å, 3 µm, 150 x 4 mm)

Eluent

70% Water
(+ 10 mM Borax + 10 mM Soda)
30% Acetonitrile

Flow


1 ml/min

Detection

UV, 230 nm

Fig. 5: RP Separation of Indole Alkaloids

pH 9,5



Column

Sachtpore RP (300 Å, 5 µm, 150 x 4 mm)

Eluent

70% Water
(+ 10 mM KH₂PO₄ 30 mM Na₂CO₃)
30% Acetonitrile

Flow


1 ml/min

Detection

UV, 280 nm

Fig. 6: RP Separation of Antihypertensives

pH 9,5



Column

Sachtpore RP (300 Å, 5 µm, 150 x 4 mm)

Eluent

70% Water
(+ 40 mM NH₃ + 10 mM K₂HPO₄
+ 10 mM Na₂CO₃)
30% Acetonitrile

Flow


1 ml/min

Detection

UV, 230 nm

Fig. 7: RP Separation of Ergot Alkaloids

pH 9,5



Column

Sachtopore RP (200 Å, 5 µm, 150 x 4 mm)

Eluent

70% Water
 (+ 10 mM K₂HPO₄ + 30 mM Na₂CO₃)
 30% Acetonitrile

Flow


1 ml/min

Detection

UV, 220 nm

Fig. 8: RP Separation of Antipyretic Drugs

pH 9



Column

Sachtopore RP (300 Å, 5 µm, 150 x 4 mm)

Eluent

91% Water
 (+ 5 mM Borax)
 9% Acetonitrile

Flow


1 ml/min

Detection

UV, 230 nm

Fig. 9: RP Separation of Antihypertensives

pH 9



Column

Sachtopore RP (300 Å, 3 µm, 150 x 4 mm)

Eluent

81% Water
 (+ 10 mM Na₂B₄O₇ + 1 mM H₃BO₃)
 19% Acetonitrile

Flow


1 ml/min

Detection

UV, 220 nm

Fig. 10: RP Separation of Antiphlogistic Drugs

pH 2,3



Column

Sachtopore RP (300 Å, 3 µm, 150 x 4 mm)

Eluent

74% Water
 (+ 50 mM H₃PO₄ + 25 mM KH₂PO₄)
 26% Acetonitrile

Flow


1 ml/min

Detection

UV, 280 nm

Fig. 11: RP Separation of Antiphlogistic Drugs

pH 2,3



Column

Sachtopore RP (300 Å, 3 µm, 150 x 4 mm)

Eluent

74% Water
(+ 50 mM H₃PO₄ + 25 mM KH₂PO₄)
26% Acetonitrile

Flow


1 ml/min

Detection

UV, 220 nm

Fig. 12: RP Separation of Antiphlogistic Drugs

pH 2,2



Column

Sachtopore RP (200 Å, 5 µm, 150 x 4 mm)

Eluent

73% Water
(+ 40 mM NaH₂PO₄ + 55 mM H₃PO₄)
27% Acetonitrile

Flow


1 ml/min

Detection

UV, 230 nm

Fig. 15: RP Separation of Antiseptics

pH 1,5



Column

Sachtopore RP (200 Å, 5 µm, 150 x 4 mm)

Eluent

85% Water
(+ 20 mM H₃PO₄)
15% Acetonitrile

Flow


0,8 ml/min

Detection

UV, 210 nm

Fig. 16: RP Separation of Antiseptics

pH 1,5



Column

Sachtopore RP (200 Å, 5 µm, 150 x 4 mm)

Eluent

75% Water
(+ 20 mM H₃PO₄)
25% Acetonitrile

Flow

1 ml/min

Detection

UV, 210 nm